

# Spore News

Volume 3, Number 2  
June 2006



Boden Larson  
IT Manager

## Z-Value Calculation

### What is a Z-value?

A Z-value is defined as the number of degrees (Celsius or Fahrenheit) required to change a D-value by one factor of ten. In the practical sense, it is a measure of how susceptible a spore population is to changes in temperature. For example, if the Z-value of a population is 10 degrees, then increasing the sterilization temperature 10 degrees will result in a log reduction of the D-value.

The Z-value can be found by plotting D-values against temperatures on a semi-logarithmic scale and adding a line that best fits the data. The absolute value of the reciprocal of the slope of this line will be the Z-value.




Figure 1, D-values at 121, 124, and 129C plotted against Temperature

Since D-values are plotted logarithmically, the line of best fit will be of exponential form and will therefore appear linear on semi-logarithmic paper, as shown above in Figure 1.

### Why use a line-of-best-fit?

In a perfect world, the D-value of a particular spore will decrease exponentially as temperature is increased. In reality however, D-values will not decrease *perfectly* due to natural variations and experimental error. A line-of-best-fit is a statistically correct method of representing the data points.

### How do I determine a Z-value graphically?

Calculating a Z-value graphically is relatively easy. Plot at least 3 D-value/temperature pairs on a semi-logarithmic graph and, using your best judgement, draw a straight line through the points that most closely fits the data. From this line, determine the number of degrees required to change the D-value by one factor of 10 (see Figure 2 below).




Figure 2, Calculating a Z-value graphically.

In the above case the Z-value appears to be approximately 8°C.

## Can I determine a Z-value mathematically?

Yes. Estimating a Z-value graphically is far less perfect than calculating a Z-value mathematically.

To accomplish this we must know the equation for our line-of-best-fit, or more specifically the slope of our line, which can be found using a log-linear regression analysis.

The following normal equations can be used to find the equation of a linear regression line:

$$\sum_{i=1}^n y_i = b \cdot n + m \cdot \sum_{i=1}^n x_i \quad (\text{eq 1})$$

$$\sum_{i=1}^n x_i y_i = b \cdot \sum_{i=1}^n x_i + m \cdot \sum_{i=1}^n x_i^2 \quad (\text{eq 2})$$

where  $b$  is the y-intercept of the regression line

$m$  is the slope

and  $n$  is the number of  $x, y$  points.

Since, however, we are dealing with  $x, \log y$  pairs, the above equations must be modified to reflect our log-linear scale:

$$\sum_{i=1}^n \log y_i = b \cdot n + m \cdot \sum_{i=1}^n x_i \quad (\text{eq 3})$$

$$\sum_{i=1}^n x_i \cdot \log y_i = b \cdot \sum_{i=1}^n x_i + m \cdot \sum_{i=1}^n x_i^2 \quad (\text{eq 4})$$

Now, using the method of least squares, we can find the following formula for the slope of our line:

$$m = \frac{n \cdot \sum_{i=1}^n (x_i \cdot \log y_i) - \sum_{i=1}^n x_i \cdot \sum_{i=1}^n \log y_i}{n \cdot \sum_{i=1}^n x_i^2 - (\sum_{i=1}^n x_i)^2} \quad (\text{eq 5})$$

where  $n$  is the number of D-value/temperature pairs (data points)  
 $x$  represents temperature  
and  $y$  represents D-value

From this equation we can find the Z-value by taking the absolute value of the reciprocal of this slope, as such:

$$Z = \left| \frac{1}{m} \right| \quad (\text{eq 6})$$

**Why is the Z-value the reciprocal of the slope? Explain further.**

Take another look at Figure 2. As you can see, the slope ("rise over run") of our regression line is:

$$m = \frac{\Delta D \text{value}}{\Delta T \text{emperature}} \Rightarrow \frac{\Delta D}{\Delta T} \quad (\text{eq 7})$$

By the definition of a Z-value, we are only concerned with a 1 log reduction in the D-value. Therefore we fix our change in D-value ( $\Delta D$ ) at 1.

$$m = \frac{1}{\Delta T} \quad (\text{eq 8})$$

In so doing, this change in temperature now becomes our Z-value.

$$m = \frac{1}{Z} \quad (\text{eq 9})$$

And solving for Z gives us:

$$Z = \frac{1}{m} \quad (\text{eq 10})$$

Finally, since the Z-value should always be a positive value (and slope of our regression line can be negative), we take the absolute value of our equation:

$$Z = \left| \frac{1}{m} \right| \quad (\text{eq 6})$$

**Example 1:**

Calculate a Z-value from the following D-values:

Temperature	D-value
115C	8.9 minutes
121C	2.1 minutes
124C	0.9 minutes

Construct a table to simplify the calculation:

D-Values (y)	Temperature (x)	log(y)	x <sup>2</sup>	x(log(y))
8.9	115	0.9494	13225	109.181
2.1	121	0.3222	14641	38.9862
0.9	124	-0.0458	15376	-5.6792
-----	Σx = 360	Σ log(y) = 1.2258	Σx <sup>2</sup> = 43242	Σx(log(y)) = 142.488

number of x, y pairs = n = 3

Now plug these values into the equation for the slope of our regression line (eq 5).

$$m = \frac{3 \cdot 142.4880 - 360 \cdot 1.2258}{3 \cdot 43242 - 360^2} = -13.824/126 = -0.1097$$

Now solve for Z (eq 6).

$$Z = \left| \frac{1}{-0.1097} \right| = 9.1 \square$$

**Ok, I understand Z-values. How do I determine a D-value from a Z-value and known D-values?**

While this task is self-explanatory when done graphically, it is rather difficult to perform mathematically.

First we must define the problem. If we have, say, three D-value/temperature pairs, then we are looking for a fourth D-value at a *given* temperature. Therefore we are solving for an unknown D at a specified temperature. The unknown D-value should fit our regression line in such a way as not to affect its slope, thereby not changing the Z-value.

Performing this task will require that we know the entire data set that was used to generate the Z-value specified. While all of this data may not be included on a certificate of performance, it can in most cases be obtained by contacting SGM Biotech, Inc.

We must expand our original equation for the slope of a log-linear regression line to include the unknown D-value.

$$m = \frac{n \cdot \left( \sum_{i=1}^{n-1} [x_i \cdot \log y_i] + x_n \cdot \log y_n \right) - \sum_{i=1}^n x_i \cdot \left( \sum_{i=1}^{n-1} [\log y_i] + \log y_n \right)}{n \cdot \sum_{i=1}^n x^2 - \left( \sum_{i=1}^n x \right)^2} \quad (\text{eq 11})$$

Since the Z-value is known, we can substitute for m:

$$\frac{-1}{Z} = \frac{n \cdot \left( \sum_{i=1}^{n-1} [x_i \cdot \log y_i] + x_n \cdot \log y_n \right) - \sum_{i=1}^n x_i \cdot \left( \sum_{i=1}^{n-1} [\log y_i] + \log y_n \right)}{n \cdot \sum_{i=1}^n x^2 - \left( \sum_{i=1}^n x \right)^2} \quad (\text{eq 12})$$

Note that we've substituted -1/Z for m instead of the absolute value of 1/Z. This is due to the fact that our D-values are assumed to decrease as temperature increases, thereby causing the slope of our regression line to be negative. If the D-value was to increase with temperature then we would use 1/Z, however it is highly unlikely that this would ever be the case.

So, what are we solving for? Well, we know Z; we know  $x_1, x_2 \dots x_n$  as our temperatures; and we know  $y_1, y_2 \dots y_{n-1}$ . What we don't know is  $y_n$ , which is of course the unknown D-value.

Solving for  $y_n$ , our unknown D-value, yields:

$$y_n = 10^{\frac{\left[ \frac{-n \sum_{i=1}^n x_i^2 + \left( \sum_{i=1}^n x_i \right)^2 \right] \cdot Z - n \sum_{i=1}^{n-1} (x_i \log y_i) + \sum_{i=1}^n x_i \cdot \sum_{i=1}^{n-1} \log y_i}{n x_n - \sum_{i=1}^n x_i}}}{Z}} \quad (\text{eq 13})$$

While this equation is somewhat cumbersome, its difficulty can be reduced using the methods of the following example.

### Example 2:

Calculate a D-value for 129C from the following data:

Temperature	D-value
115C	8.9 minutes
121C	2.1 minutes
124C	0.9 minutes

$$Z\text{-value} = 9.5C$$

Construct a table to simplify the calculation:

D-Values (y)	Temperature (x)	log(y)	x <sup>2</sup>	x(log(y))
8.9	115	0.9494	13225	109.181
2.1	121	0.3222	14641	38.9862
0.9	124	-0.0458	15376	-5.6792
y <sub>n</sub>	129	-----	16641	-----
-----	Sum of x <sub>1</sub> ..x <sub>n</sub> = 489	Sum of log(y)...log(y <sub>n-1</sub> ) = 1.2258	Sum of x <sub>1</sub> <sup>2</sup> ..x <sub>n</sub> <sup>2</sup> = 59883	Sum of x <sub>1</sub> (log(y <sub>1</sub> ))..x <sub>n-1</sub> (log(y <sub>n-1</sub> )) = 142.488

number of x values = n = 4

$$x_n = x_4 = 129$$

Now plug these values into the equation (eq 13).

$$y_4 = 10^{\frac{\frac{-4.59883 + 239121}{9.5} - 4.142.488 + 489 \cdot 1.2258}{4 \cdot 129 - 489}} = 0.31 \text{ minutes } \square$$

We can check this value by including it in our data set and calculating a new Z-value. The new Z-value should be the same as the one stated in the problem, 9.5C.

### Can I calculate a D-value using only the data on the certificate?

In most cases a certificate of performance will at least include the D-value at 121C and the Z-value (where applicable). While it is possible to estimate a D-value from this data, it is not recommended.

Since the Z-value only tells us about the slope of the original regression line and not its y-intercept, the line can exist anywhere along the y-axis. While it will always generate the same Z-value, the D-values it reports will be skewed.

Because we don't know where exactly to place our line, we're simply going to make the assumption that it passes through the D-value at 121C. Remember that in a perfect world

all D-values would exist on our regression line, as they would decrease perfectly exponentially.

For our solution, we're going to use the most basic form for the slope of our line:

$$m = \frac{y_2 - y_1}{x_2 - x_1} \quad (\text{eq 14})$$

Which in this case will be:

$$m = \frac{\log D_{121} - \log D_x}{121 - T_x} \quad (\text{eq 15})$$

where  $D_{121}$  is the D-value at 121C

$D_x$  is the unknown D-value

and  $T_x$  is the temperature at which the unknown D-value exists.

Solving for  $D_x$  results in:

$$D_x = 10^{\frac{1}{Z}(121 - T_x) + \log D_{121}} \quad (\text{eq 16})$$

IFF  $T_x < 121$

$$D_x = 10^{\frac{-1}{Z}(T_x - 121) + \log D_{121}} \quad (\text{eq 17})$$

IFF  $T_x > 121$ .

### What do the standards say about Z-values?

AAMI, ISO, and EN all discuss Z-values.

	EN 866-3	ISO 11138-3	AAMI ST34	SGM Biotech
Valid Temperature Range	110-130C	110-130C	Not stated.	<b>110-130C</b>
Number of D-values Required	2	2	2	<b>3</b>
Minimum Z-value	6C	6C	Not stated.	<b>6C</b>
Increment (accuracy)	Not stated.	0.1C	Not stated.	<b>0.1C</b>
Formula	Not stated.	$Z = \frac{T_2 - T_1}{\log D_1 - \log D_2}$	Graphically. Regression analysis acknowledged.	<b>Log-linear regression. See equation 6.</b>



## **What about the USP?**

The United States Pharmacopoeia does not discuss Z-values as a performance statistic.

## **Why are SGM Biotech's standards more strict than EN, ISO, and AAMI?**

As every scientist and statistician knows, the more data used in an estimation, the more accurate that estimation will be. We do not feel that using two D-values to determine a Z-value is proper, since the slope of a regression line through only two points will vary greatly from a line through three or more points. Three D-values is the minimum number of points required for practical yet accurate Z-value estimates.

## **I can't find the Z-value on the certificate that came with the product!?**

Z-values are generally only reported for products used in steam and dry heat processes.

**Please email us with topics you would like to see addressed in “Spore News”.**

