

Prepared Culture Media

Soybean Casein Digest Broth is the recommended media (USP 24) for use with Biological Indicators. Use of another media may detrimentally affect outgrowth of spore. Mesa's TSB is manufactured and tested to meet USP and NCCLS sterility and growth promotion guidelines. It is also specially formulated for use with our spore products.

Prepared Media is packaged in Type I borosilicate glass tubes with polypropylene screw-caps (white rubber liners). Available sizes: 13mm X 100mm, 16mm X 100mm, and 20mm X 100mm test tubes. Also available in 500ml and 1L screw-capped, pourable bottles. Other sizes available upon request.



24-Hour Incubation

A reduced incubation time of 24 hours at 60°C +/- 2°C has been validated for steam sterilization when the Mesa Bacterial Spore Strip (either the single species *G. stearothermophilus* spore strip or the dual species spore strip) is used in conjunction with Mesa's (OMF) Modified Tryptic Soy Broth with Bromocresol Purple. 510k information available upon request.



Dehydrated Culture Media

Mesa recommends Acumedia TSA and TSB for use in verifying the certified label claims for our products. Acumedia has an extensive line of high quality dehydrated culture media at affordable prices. Mesa uses Acumedia in-house for the manufacture and QC testing of our products. The entire Acumedia line is now available through Mesa Labs and many of our international distributors.

Ordering Information

Tryptic Soy Broth	Catalog Number	Fill	Quantity	Tube Size
Tryptic Soy Broth	TSB-C16	6.5 ml	100 / box	16mm x 100mm
Tryptic Soy Broth w/Phenol Red	TSB-PR16	6.5 ml	100 / box	16mm x 100mm
Tryptic Soy Broth w/Bromocresol Purple	TSB-BP16	6.5 ml	100 / box	16mm x 100mm
Tryptic Soy Broth in Screw Capped Bottle	TSB-1000	1000 ml	minimum purchase 10 bottles	