

Mesa Biological Indicators

Spore Strips and Culture Media

MagnaAmp®

MagnaAmp® is a self-contained biological indicator ampoule for monitoring steam sterilization of liquids in containers or washer sterilizers. The MagnaAmp® biological indicator ampoules contains *Geobacillus stearothermophilus* spores and a specially-formulated culture medium that turns a dramatic yellow if spores grow.

The MagnaAmp® provides easy visual confirmation of sterilization within 48-hours, without the need for difficult and time-consuming laboratory work.



- 48-hour results
- Small (1ml), conveniently shaped ampoule
- Designed to sink in container
- Negative controls are included with each box of product
- No laboratory required
- Dramatic color change

Storage

Refrigerate at 2-8°C.

Certification

MagnaAmp® biological indicators are traceable to a recognized culture collection and certified for population, D-value, z-value, F₀ value and performance.

Ordering Information

Species	Catalog #	Min. Population	Sets/Box
<i>G. stearothermophilus</i>	MA/5	1.0 x 10 ⁵	50
<i>G. stearothermophilus</i>	MA/6	1.0 x 10 ⁶	50