

SterilAmp® II is a compact, self-contained biological indicator for monitoring industrial steam sterilization of liquids in small containers. SterilAmp® II precisely integrates the process lethality.

- **48-hour results**
- Small, conveniently-shaped ampoule
- No laboratory required
- Dramatic color change
- Packaged with 10 negative controls

26mm x 6.5mm



PASS

FAIL

SterilAmp® II contains a specially formulated culture medium that turns a dramatic yellow when spores grow and provides easy visual confirmation of sterilization within 48 hours, without the need for difficult and time-consuming laboratory work.

Storage

Refrigerate at 2-8°C.

Certification

SterilAmp® biological indicators are traceable to a recognized culture collection and certified for population, D-value, z-value, Fo value, and performance.

Ordering Information

Species	Catalog #	Min. Population	Qty/Box
<i>G. stearothermophilus</i>	SA/4	1.0 x 10 ⁴	100 Tests / 10 NC
<i>G. stearothermophilus</i>	SA/5	1.0 x 10 ⁵	100 Tests / 10 NC
<i>G. stearothermophilus</i>	SA/6	1.0 x 10 ⁶	100 Tests / 10 NC