

Apex Product List

Description	Reorder #	Quantity/ Pkg
Tyvek/Tyvek Packaged Biological Indicator, <i>G. stearothermophilus</i> (#12980) inoculated stainless steel discs, 10 ⁶ / carrier.	HMV-091	100
Tyvek/Tyvek Packaged Biological Indicator, <i>G. stearothermophilus</i> (#12980) inoculated stainless steel discs, 10 ⁶ / carrier - packaged in triplicate.	HMV-091T	33 pkgs of 3 (99 discs)
Tyvek/Tyvek Packaged Biological Indicator, <i>G. stearothermophilus</i> (#12980) inoculated stainless steel discs, 10 ⁴ / carrier.	KCD-404	100
Bare Metal Biological Indicator, <i>G. stearothermophilus</i> (#12980) inoculated stainless steel ribbon, 10 ⁶ / carrier.	SBC-327	100
Tri-Scale Biological Indicator -- Tyvek packaged BI with 10 ⁶ , 10 ⁵ , and 10 ⁴ inoculated stainless steel carriers (<i>G. stearothermophilus</i> #12980).	LOG-456	50
Tyvek/Tyvek Packaged Biological Indicator, <i>B. atrophaeus</i> (#9372) inoculated stainless steel discs, 10 ⁴ -10 ⁷ / carrier.	GRS series	100
Tyvek/Tyvek Packaged Biological Indicator, <i>G. stearothermophilus</i> (#7953) inoculated stainless steel discs, 10 ⁶ / carrier.	NAS-152	100
Sterile Tyvek Transfer Pockets for SBC-327.	BMWK-750	100
<i>G. stearothermophilus</i> (#12980) suspension calibrated for gaseous hydrogen peroxide, (10 ⁷ / 0.1 ml).	LAU-196	10 ml vial

Contact us directly for: Custom Manufactured Biological Indicators, Spore Crop Production from Environmental Isolates, Resistance Testing (Steam, H₂O₂) and Sterilization Process Investigations.

For more information, contact Mesa Labs or your local distributor

Ordering Information:

Tel: 303-987-8000

Fax: 406-585-9219

apex-cust-service@mesalabs.com

10 Evergreen Drive
Bozeman, MT 59715 USA
www.mesalabs.com

Payment Terms: Net 30 days WAC, MasterCard, Visa, AmEx accepted

Shipping: Prepaid and added to invoice

Product specifications and pricing subject to change without notice.