


The Mesa 1510 incubator is designed for convenient and effective incubation of Mesa self-contained biological indicators at 57°C including MagnaAmp® and ProSpore Ampoule 1mL liquid submersible biological indicators.



1510 SCBI Incubator

Convenient

- Simultaneous incubation of up to 15 BIs in numbered incubation cavities.
- Multinational 110V/220V power adapter with interchangeable blades.
- Digital Temperature Display

User Friendly

- Small footprint.
- Fast warm up time (30-45 minutes).
- Designed for 24-hour-a-day use.
- Reliable operation.

Certifications

RoHS compliant, CSA C/US and CE Mark