

The EZTest® Biological Test Pack is a small, convenient, FDA 510(k) notified product containing EZTest® self-contained biological indicators. Designed to parallel the performance of the AAMI Test Pack, each pack includes a record card with a chemical integrator.

- FDA 510(k) notified (steam and EO)
- 24-hour results / 10-hour results (SmartRead EZTest)
- Easy to culture
- No laboratory required
- Dramatic color change
- Release loads immediately*



Storage

Store at room temperature.

Certification

EZTest® Biological Test Packs are traceable to a recognized culture collection and certified for population, and D-value, z-value (steam) and performance.

Ordering Information

Steam Biological Indicator Test Pack (PCD)			
Species	Catalog #	Min. Population	Qty/Box
<i>G. stearothermophilus</i>	STP/25	1.0 x 10 ⁵	25 Tests / 25 Controls
Smart-Read® EZTest® Steam Biological Indicator Test Pack (PCD)			
Species	Catalog #	Min. Population	Qty/Box
<i>G. stearothermophilus</i>	SRSP/05	1.0 x 10 ⁵	25 Tests / 5 Controls
<i>G. stearothermophilus</i>	SRSP/25	1.0 x 10 ⁵	25 Tests / 25 Controls

**Non-implant loads only, loads containing implants should be quarantined until biological results are known. Instant readout card may be used to aid in the release of implant loads prior to biological results in emergency situations only.*