

Mesa Biological Indicators

1410 Adjustable SCBI Incubator

1410 Adjustable Self-Contained Biological Indicator (SCBI) Incubator

The Mesa 1410 incubator is designed for convenient and effective incubation of Mesa self-contained biological indicators at 37°C (*Bacillus atrophaeus*) or 57°C and 60°C (*Geobacillus stearothermophilus*).



The incubator features a built-in thermometer with LED display, 13 incubation cavities and an integrated activation (crushing) cavity.

Convenient

- Simultaneous incubation of up to 13 BIs in numbered incubation cavities.
- Multinational 110V/220V power adapter with interchangeable blades.
- Incubation temperature adjustable for BIs containing *B. atrophaeus* (37°C) or *G. stearothermophilus* (57°C and 60°C).
- Supports 10-hour Smart-Read[®] BIs at 60°C.

User Friendly

- Integrated activation (crushing) cavity.
- Very fast warm up time (10-15 minutes depending on selected temperature).
- Designed for 24-hours-a-day use.
- Integrated thermometer with digital display.

Certifications

RoHS compliant, CSA C/US and CE Mark.