

Biological and Chemical Indicator

Results Log Sheet

Date	Time	Processor Serial Number and/or Cycle #	Location and/or Load	Operator	Incubator Temp	Lot Number SS (spore strips) CM(culture media)	Results		Notes	Date/Time Out Initials			
							BI	CI					
							Control	Test		✓	✗		
								24 hours	Final				
						S.S.	+ -	+ -	+ -				
						C.M.							
						S.S.	+ -	+ -	+ -				
						C.M.							
						S.S.	+ -	+ -	+ -				
						C.M.							
						S.S.	+ -	+ -	+ -				
						C.M.							
						S.S.	+ -	+ -	+ -				
						C.M.							
						S.S.	+ -	+ -	+ -				
						C.M.							
						S.S.	+ -	+ -	+ -				
						C.M.							
						S.S.	+ -	+ -	+ -				
						C.M.							
						S.S.	+ -	+ -	+ -				
						C.M.							
						S.S.	+ -	+ -	+ -				
						C.M.							
						S.S.	+ -	+ -	+ -				
						C.M.							

Biological

- + Turbidity and/or color change to or towards yellow
- No turbidity and no color change (sterilization conditions achieved)

Chemical

- ✓ Color Loss Complete (sterilant potency correct)
- ✗ Color Loss Incomplete