

ProChem Chemical Indicators

Emulators, Integrators and Process Indicators

From indicator labels to Class 6 emulators, Mesa has a full line of chemical indicators for healthcare and industrial users alike. Mesa offers indicators for a variety of processes including Steam, EO, Gamma, Plasma and Dry Heat.

Integrators and Emulators

ProChem SSW - Steam Class 5

- Moving front Integrator, no interpretation necessary.
- Effective across the normal sterilization range.
- Parallels the performance of a biological indicator.




Catalog # CI-SSW
250 Indicators / box
3 year expiry

ProChem EXT - Steam Class 5

- 7 minutes at 135 °C / 20 minutes at 121 °C
- Laminated
- Color drives darker than color standard
- Permanent color change

Catalog # CI-EXT
250 Indicators / box
5 year expiry



Process Indicators

ProChem SSI - Steam Class 4

International Use Only

- 3 minutes at 132 °C - 135 °C
- Most cost effective CI solution
- Drives darker than the color standard
- Indicator Perforated; Indicator spot / LOT# printed on both sides

Catalog # CI-SINT
250 Indicators / box
3 year expiry



ProChem OK - Steam Sterilization Multi-Spot Indicator

International Use Only

Catalog # CI-OK121
250 Indicators / bag
3 year expiry

- Pre-vacuum sterilizers operating at 121 °C and 134 °C
- Multiple spots with varied timing
- Single indicator suitable for multiple cycle lengths



ProChem Dry Heat Indicator Labels

- Green to dark brown / black color change
- Distinguishes between processed and unprocessed items


Catalog # CI-DHI
1,000 indicators / roll
3 year expiry

Bowie-Dick Type Tests


- Distinct, non-reversible pink to black color change
- Meets AAMI, EN 867-3 and ISO 11140 air removal
- Preassembled, disposable test pack

ProChem B-D

Catalog # CI-BDTLF-20
Catalog # CI-BDTLF




Box of 20
Case of 80



ProChem B-D 121*

Catalog # CI-BDTLF-121-20
Catalog # CI-BDTLF-121




ProChem BD110*

Individual Test Sheets

- Distinct, non-reversible blue to black color change
- Meets AAMI, EN 867-3 and ISO 11140 air removal

Catalog # CI-BD110
100 sheets / package
3 year expiry



**Industrial/International Use Only*