

From indicator labels to Class 6 emulators, Mesa has a full line of chemical indicators for healthcare and industrial users alike. Mesa offers indicators for a variety of processes including Steam, EO, Gamma, Plasma and Dry Heat.

Integrators and Emulators

ProChem SSW - Steam Class 5

- Moving front Integrator, no interpretation necessary.
- Effective across the normal sterilization range.
- Parallels the performance of a biological indicator.




Catalog # CI-SSW
250 Indicators / box
3 year expiry

ProChem EXT - Steam Class 5

- 7 minutes at 135°C / 20 minutes at 121°C
- Laminated
- Color drives darker than color standard
- Permanent color change


Catalog # CI-EXT
250 Indicators / box
5 year expiry



ProChem SSE - Steam Class 6

- 3 minutes at 135°C / 4 minutes at 132°C
- Printed color standard
- Laminated
- Lead Free

Catalog # CI-SSIX
250 Indicators / bag
18 month expiry



Process Indicators

ProChem SSI - Steam Class 4

International Use Only

- 3 minutes at 132 - 135°C
- Most cost effective CI solution
- Drives darker than the color standard
- Indicator Perforated; Indicator spot / LOT# printed on both sides

Catalog # CI-SINT
250 Indicators / box
3 year expiry



ProChem OK - Steam Sterilization Multi-Spot Indicator

International Use Only



- Pre-vacuum sterilizers operating at 121°C and 134°C
- Multiple spots with varied timing
- Single indicator suitable for multiple cycle lengths

Catalog # CI-OK121
250 Indicators / bag
3 year expiry



ProChem Dry Heat Indicator Labels

- Green to dark brown / black color change
- Distinguishes between processed and unprocessed items

Catalog # CI-DHI
1,000 indicators / roll
3 year expiry

Gamma Indicator Labels

- Yellow to red color change
- 2 year expiry



Catalog # CI-R02 (METO® gun label)
1,000 indicators / roll
Catalog # CI-R03 (12.7mm dots)
5,000 indicators / roll

Ethylene Oxide Indicator Labels

- Purple to green color change
- 2 year expiry



Catalog # CI-E02 (METO® gun label)
1,500 indicators / roll
Catalog # CI-E03 (12.7mm dots)
5,000 indicators / roll

Bowie-Dick Type Tests

- Distinct, non-reversible pink to black color change
- Meets AAMI, EN 867-3 and ISO 11140 air removal
- Preassembled, disposable test pack




ProChem B-D

Catalog # CI-BDTLF-20
Catalog # CI-BDTLF

Box of 20
Case of 80

ProChem B-D 121*

Catalog # CI-BDTLF-121-20
Catalog # CI-BDTLF-121




ProChem BD110*

Individual Test Sheets

- Distinct, non-reversible blue to black color change
- Meets AAMI, EN 867-3 and ISO 11140 air removal

Catalog # CI-BD110
100 sheets / package
3 year expiry



***Industrial/International Use Only**