

EZTest® is a self-contained biological indicator for monitoring steam, ethylene oxide or hydrogen peroxide sterilization. EZTest contains bacterial spores on a small carrier packaged within a small, thermoplastic culture tube. Inside the thermoplastic culture tube is a sealed-glass ampoule of specially-formulated soybean casein digest culture medium containing a color indicator which turns a dramatic yellow when spores grow. EZTest is easy-to-use and produces visual results within 24 hours (steam or hydrogen peroxide) or 48 hours (ethylene oxide) without laboratory transfers.



- FDA 510(k) notified (steam and EO)
- 24-hour results for steam or H2O2
- 48-hour results for EO
- Easy to culture
- No laboratory required
- Dramatic color change

Ordering Information

| EZTest® for Steam / Flash | | | |
|--------------------------------------|-----------|-----------------------|---------|
| Species | Catalog # | Min. Population | Qty/Box |
| <i>G. stearothermophilus</i> | EZS/5 | 1.0 x 10 ⁵ | 100 |
| <i>G. stearothermophilus</i> | EZS/525 | 1.0 x 10 ⁵ | 25 |
| <i>G. stearothermophilus</i> | EZS/6 | 1.0 x 10 ⁶ | 100 |
| <i>G. stearothermophilus</i> | EZS/625 | 1.0 x 10 ⁶ | 25 |
| EZTest® for Ethylene Oxide | | | |
| Species | Catalog # | Min. Population | Qty/Box |
| <i>B. atrophaeus</i> | EZG/6 | 1.0 x 10 ⁶ | 100 |
| <i>B. atrophaeus</i> | EZG/625 | 1.0 x 10 ⁶ | 25 |
| EZTest® for Hydrogen Peroxide (H2O2) | | | |
| Species | Catalog # | Min. Population | Qty/Box |
| <i>G. stearothermophilus</i> | EZH/5 | 1.0 x 10 ⁵ | 100 |
| <i>G. stearothermophilus</i> | EZH/6 | 1.0 x 10 ⁶ | 100 |