

Apex BI Product List

Mesa Labs Vaporized Hydrogen Peroxide Biological Indicators

APEX Biological Indicator products are designed specifically for the pharmaceutical, food and medical device industries utilizing H₂O₂ vapor sterilization. The stainless steel carrier is designed for use with today's isolator and filling line applications.

| Item | Reorder Number | Population | ATCC Number | Package Configuration |
|--|----------------|---------------------|-------------|-----------------------|
| Stainless Steel Discs in Tyvek/Tyvek <i>G. stearothermophilus</i> | HMV-091E3 | 10x3 | 12980 | 100/pkg |
| | HMV-091 | 10x6 | 12980 | 100/pkg |
| | KCD-404 | 10x4 | 12980 | 100/pkg |
| Bare Metal BI Stainless Steel Ribbon <i>G. stearothermophilus</i> | SBC-327 | 10x6 | 12980 | 100/pkg |
| Tri-Scale BI; Tyvek pkg BI w/ 10x6, 10x5, 10x4 Stainless Steel Carriers <i>G. stearothermophilus</i> | LOG-456 | 10x4, 10x5, 10x6 | 12980 | 50/pkg |
| Stainless Steel Discs in Tyvek/Tyvek <i>B. atrophaeus</i> | GRS-090E3 | 10x3 | 9372 | 100/pkg |
| | GRS-090 | 10x6 | 9372 | 100/pkg |
| Stainless Steel Discs in Tyvek/Tyvek <i>G. stearothermophilus</i> | NAS-152E3 | 10x3 | 7953 | 100/pkg |
| | NAS-152 | 10x6 | 7953 | 100/pkg |
| Sterile Tyvek Transfer Pockets for SBC-327 | BMWK-750 | | | 100/pkg |
| Spore Suspension for gaseous hydrogen peroxide <i>G. stearothermophilus</i> | LAU-196 | 10x7/0.1 ml | 12980 | 10 ml vial |
| | LPT-606 | 10x7/0.1 ml | 7953 | 10 ml vial |
| Releasat Purple, for use with HMV-091, KCD-404 and LOG-456 3.6 Fill Volume in 16mm x 100mm Tubes | PM/100 | -- | -- | 100 units |

For more information, contact Mesa Labs or your local distributor.

