

SUPPLY INTERRUPTION NOTICE

February 10, 2020

Notification ID: CN-200201

Notification Type: Supply Interruption – PCD 5.x Products

Dear Valued Customer,

This letter is to inform you of a supply interruption from Mesa's 625 Zoot Way Bozeman, MT Manufacturing Facility for the following products:

- PCD 5.1, PCD 5.2, PCD 5.5, PCD 5.13, PCD 5.62, PCD 5.92

This interruption of supply is a result of a material shortage due to the newest receipt of Type 5 pouches failing to meet Mesa's quality control standards. At this time, we do not have a projected date as to when the manufacture of these products will resume. We are working with our supplier to obtain lead times to produce the raw film and conversion of the film to pouches. The table below contains information for alternative products that have comparable resistance values:

Current Products		Alternative Products			
PCD Type	Nominal D-value	PCD Type	Nominal D-value	PCD Type	Nominal D-value
PCD 5.1	4minutes	PCD 3.5	4minutes	PCD 3.6	5minutes
PCD 5.2	4minutes	PCD 3.5	4minutes	PCD 3.6	5minutes
PCD 5.5	6minutes	PCD 3.6	5minutes	PCD 7.1 PCD 7.2	9minutes
PCD 5.13	13minutes	PCD 6.1 PCD 6.2	12minutes	PCD 7.5 PCD 7.7	14minutes
PCD 5.62	13minutes	PCD 6.1 PCD 6.2	12minutes	PCD 7.5 PCD 7.7	14minutes
PCD 5.92	13minutes	PCD 6.1 PCD 6.2	12minutes	PCD 7.5 PCD 7.7	14minutes


As soon as we have dates for the manufacture of the film and conversion to pouches or if the status of this situation changes, we will update this information. We appreciate your business and apologize for the inconvenience this condition has created for you.

Please contact us with any questions or concerns at bi-support@mesalabs.com



Robert Bradley – Sr. Director – Bozeman Operations

10Feb2020
Date



Nicole DeWald – Manager – Quality Assurance

10Feb2020
Date